


# MEASUREMENT GUIDE

It really is very easy to take measurements for your new stairlift. All you will need is a good quality metal tape measure. First, take a general overview of your stairway layout and decide which side of the staircase the stairlift will be mounted on. Be aware that the stairlift rail will overhang the top landing by 7" (20cm) and meet the floor at the foot of the stairs 13" (33cm) beyond the 'nose' of the first step. Take note of any fixtures, such as window sills, hand rails, central heating radiators, doorways etc. that intrude into the stairway area - or could be blocked by the presence of the stairlift. Next, simply stand at the top of the staircase and feed the tape measure down the steps until it goes out at an angle and hits the floor. What is the measurement at your toes?

## Step 1 - Staircase length

Please provide the following measurement:



## ADDITIONAL MEASUREMENTS NEEDED

- 1) If there is a door at the top – what is the distance from the face of the door to the edge of the landing when the door is closed?
- 2) Standing at the bottom and looking up, which side of the staircase will your stairlift be mounted on?
- 3) Is there 24 inches of clearance at the bottom of the staircase to allow the chair to fully lower to the floor?